

# ACCESSORIES

The right accessories make it easier to work with autoclaves and optimise the loading volume as well as the cycle times. Baskets in different versions maximise loading to suit the chamber size of the respective unit. Loading the autoclave with heavy sterilizing goods is significantly facilitated by lifting devices or loading systems.

## Accessories for LABOKLAV free-standing units

- ◆ Baskets in various heights
- ◆ Baskets with closed bottom
- ◆ Buckets, also with lid, for waste sterilization
- ◆ Integrated magnetic stirrer
- ◆ Lifting device for heavy loads



Lifting device

## Accessories for LABOKLAV tabletop units

- ◆ Baskets in various sizes
- ◆ Trays for waste sterilization
- ◆ Middle floor mounts
- ◆ Stands
- ◆ Integrated magnetic stirrer
- ◆ Loading system

Magnetic stirrer for tabletop units



Loading system



## Accessories for all LABOKLAV laboratory autoclaves



- ◆ Condensate collection tank
- ◆ Compressors
- ◆ Cartridges for the production of demineralized water
- ◆ Integrated measurement of the quality of the demineralized water (conductivity measurement)